

DNA extraction from Oragene® Dx OGD-600 Saliva collection tubes using Isopure Buccal DNA kit

Purpose

Saliva contains buccal epithelial cells and white blood cells that contain DNA. Due to the easy access and non-invasive character, saliva is an excellent alternative to blood collection. Here we present a high-quality DNA extraction method from saliva using Isopure Buccal DNA kit. This method is applicable for extraction of DNA from Oragene Dx OGD-600 Saliva collection tubes.

Protocol

1. Sample Preparation

- a. Collect and store saliva sample according to the manufacturer's instructions

2. Lysis

- a. Incubate the Oregene tubes for 1 hour at 50°C
- b. Transfer 500 µL of sample to 2 mL 96-well plate

3. Bind

- a. Add 200 µL of Bind 1 (conditioning) to the plate
- b. Mix by pipetting up and down 10 times, or until thoroughly mixed
- c. Vortex the bottle of Bind2 to fully resuspend the beads
- d. Add 340 µL of Bind2 to the plate
- e. Mix by pipetting up and down 10 times, or until thoroughly mixed
- f. Incubate the plate for 1 minute at room temperature
- g. Place the plate on a magnet for 8 minutes (or until the supernatant is clear)
- h. Remove and discard the supernatant without disrupting the beads
- i. Remove the plate from the magnet

4. Wash

- a. Add 700 µL of 70% Ethanol Wash Buffer to plate
- b. Mix by pipetting up and down 20 times, or until thoroughly mixed
- c. Place the plate on a magnet for 2 minutes (or until the supernatant is clear)
- d. Remove and discard the supernatant without disrupting the beads
- e. Remove the plate from the magnet
- f. Repeat steps 4.a-3.e for a total of 3 washes

5. Elute

- a. Add 50 µL of EB to plate
- b. Mix by pipetting up and down 10 times, or until thoroughly mixed
- c. Place the plate on a magnet for 3 minutes (or until the supernatant is clear)
- d. Remove and Save the supernatant without disrupting the beads