



PCRCLEAN PB™

DNA Template Cleanup

Catalog Numbers: **C-1006P-5, C-1006P-50, C-1006P-250, C-1006P-450**

Contents

Protocol Manual *Revision 2.10*

Introduction	2
Process Overview	2
Product Specifications	3
Procedure - 96 Well Format	5
Procedure - 384 Well Format	7

INTRODUCTION

The PCRCLEAN PB™ purification system utilizes a proprietary magnetic bead based technology for quick and efficient, high-throughput purification of DNA from PCR and other enzyme reactions for Pac Bio instrument. PCRCLEAN PB™ utilizes an optimized buffer for selective binding of DNA, 120bp and larger to magnetic beads. The protocol mainly consists of binding, washing and elution steps. Primers, nucleotides, salts and enzymes in reaction mixture are removed during the binding and washing steps. The purified DNA product is essentially free of contaminants and suitable for Sanger sequencing or Next Gen sequencing library construction. Furthermore, primer dimers and adaptors can be effectively removed by using PCRCLEAN PB™.

PCRCLEAN PB™ can be used in the following applications:

- PCR
- Sequencing (Sanger and Next Gen library construction)
- Genotyping and mutation detection
- All enzymatic reaction clean up
- Microarray sample preparation
- Primer Walking
- Cloning and all other molecular engineering

The PCRCLEAN PB™ is adaptable to various commercially available automation platforms.

PROCESS OVERVIEW

1. Add 18ul PCRCLEAN PB™ per 10ul of PCR for DNA binding to magnetic beads.
2. Separate beads and DNA from unbound contaminants.
3. Wash beads twice with 80% Ethanol to remove salt and contaminants (Fresh 70% ethanol may be used instead).
4. Elute purified DNA from magnetic beads to a new plate.



PRODUCT SPECIFICATIONS

The PCRCLEAN PB™ kit can be performed in a tube, 96 well and 384 well formats. The following table illustrates the number of clean up reactions PCRCLEAN PB™ kits can perform depending on the sample volume of PCR.

PCRCLEAN PB™	Catalog Number
PCRCLEAN PB™ - 5 mL	C-1006P-5
PCRCLEAN PB™ - 50 mL	C-1006P-50
PCRCLEAN PB™ - 250 mL	C-1006P-250

Number of reactions:

Sample Volume (uL)	C-1006P-5 (Numbers of reactions)	C-1006P-50 (Numbers of reactions)	C-1006P-200 (Numbers of reactions)
10	277	2777	11111
20	138	1388	5555
50	55	555	2222

MATERIALS

Supplied in the Kit:

PCRCLEAN PB™ magnetic bead solution

Store at 4°C upon arrival (DO NOT FREEZE), for up to 12 months

Mix the reagent well at room temperature to completely resuspend beads prior to use. It should appear visually homogenous.

To be supplied by the User:

Apparatus

Name	Recommended Model	Recommended Vendor and P/N
96-well PCR reaction plate	96-well round-bottom microtiter plate	Corning, Inc., # 3797, www.corning.com Fisher Scientific # 07-200-105, www.fishersci.com
	96 well cycling plate	ABgene Limited, # AB-0800, AB-1000, AB-1400, www.abgene.com worldwide and Fisher Scientific www.fishersci.com in the U.S.
384-well PCR reaction plate	384 well cycling plate	Bio-Rad Laboratories, # HSP-3801, www.bio-rad.com
		ABgene product # AB-1111, www.abgene.com worldwide and Fisher Scientific www.fishersci.com in the U.S.
Magnetic PCR plate	96 well ring stand	Ambion Inc., (acquired by Applied Biosystems), # AM10050, www.appliedbiosystems.com
	384 well [MagnaBot® 384 Magnetic Separation Device]	Promega Corporation, # V8241, www.promega.com
PCR Plate Seals	Easy Peel Heat Sealing Foil	ABgene Limited, # AB-3739 and AB-3739/s, www.abgene.com worldwide and Fisher Scientific www.fishersci.com in the U.S.
Liquid handling robotics		
multichannel hand pipette		

Reagents

Reagents	Application
80% ethanol (or freshly made 70% ethanol)	Washing
10 mM TRIS-Acetate, pH 8.0	DNA elution
Reagent grade water	
10 mM Tris-Acetate pH 8.0, 0.1 mM EDTA	

PROCEDURE – 96 WELL FORMAT

1. Determine whether or not it is necessary for a plate transfer.

Note: If the PCR volume multiplied by 2.8 exceeds the volume of the PCR plate, a transfer of PCR sample to a 300 uL round bottom plate is required.

2. Gently shake the PCRCLEAN PB™ bottle to ensure a COMPLETE resuspension of any magnetic particles that may have settled. Add PCRCLEAN PB™ according to the table below:

PCR Volume (μL)	PCRCLEAN PB™ Volume at 1.8X (μL)
10	18
20	36
50	90

Note: Use the following equation to determine the volume of PCRCLEAN PB™:

$$(\text{Volume of PCRCLEAN PB™ per reaction}) = 1.8 \times (\text{PCR Volume})$$

3. Mix PCRCLEAN PB™ and PCR sample thoroughly by pipetting up and down at least 5 times.

Note: After mixing, the color of the mixture should appear homogenous.

4. Incubate the mixed samples for 5 minutes at room temperature for maximum recovery.

Note: This step allows the binding of PCR products 120bp and greater to the magnetic beads.

5. Place the reaction plate onto a 96 well magnet plate for 3 minutes or wait until the

solution is clear.

Note: Wait until the solution is clear before proceeding to the next step. Otherwise unsettled beads may become lost.

6. Aspirate the cleared solution from the reaction plate and discard.

Note: This step must be performed while the reaction plate is placed on the magnet plate. Avoid disturbing the settled magnetic beads. If beads are drawn into tips, eject for the beads to resettle and leave behind a few microliters of solution.

7. Dispense 200 μ L of 80% ethanol to each well of the reaction plate and incubate for 30 seconds at room temperature. Aspirate the ethanol and discard. Repeat for a total of two washes.

Note: It is important to perform these steps with the reaction plate on a magnetic plate. Do not disturb the settled magnetic beads. Remove all of the ethanol from the bottom of the well to avoid ethanol carryover.

A 5 min air dry at room temperature is recommended for the evaporation of the remaining traces of ethanol. **DO NOT OVERDRY THE BEADS** (bead ring appears cracked) as this will significantly decrease elution efficiency.

8. Take off the plate from the magnetic plate, add 40 μ L of elution buffer (Reagent grade water, 10mM Tris-Acetate pH 8.0, or TE buffer) to each well of the reaction plate and pipette mix 5 times.

Note: More than 40 μ L of elution buffer can be used, but using less than 40 μ L will require extra mixing to ensure the liquid comes into contact with the beads and may not be sufficient to elute the entire PCR product. Elution is quite rapid and it is not necessary for the beads to go back into solution for it to occur.

9. Place the reaction plate onto a 96 well Magnetic Plate for 1 minute to separate beads from the solution.
10. Transfer the eluate to a new plate for storage and analysis.

PROCEDURE - 384 WELL FORMAT

1. Gently shake the PCRCLEAN PB™ bottle to resuspend any magnetic particles that may have settled. Add PCRCLEAN PB™ according to the following PCR volume table:

PCR Volume (μL)	PCRCLEAN PB™ Volume (μL)
5	9
7	12.6
10	18
14	25

Note: Use the following equation to determine the volume of PCRCLEAN PB™ for a given reaction: (Volume of PCRCLEAN PB™ per reaction) = 1.8 x (PCR Volume)

Due to the constraint of the total volume of PCR plus reagent, it is not possible to purify PCR reactions larger than 14 μL with 384 well plates (14 μL reaction + 25 μL PCRCLEAN PB™ = 39 μL).

2. Mix PCRCLEAN PB™ and PCR sample thoroughly by pipette mixing 5 times and incubate at room temperature for 5 minutes.

Note: The color of the mixture should appear homogenous after mixing.

3. Place the reaction plate onto a 384 magnetic plate for 1 minute to separate the beads from solution.

Note: Wait until the solution is clear before proceeding to the next step.

4. Aspirate the cleared supernatant from the reaction plate and discard.

Note: This step should be performed while the reaction plate is placed on the 384 well magnetic plate. Do not touch the magnetic beads, which have formed a spot on the side of the well.

5. Dispense 30 µL of 80% ethanol wash solution to each well of the reaction plate and incubate for 30 seconds at room temperature. Aspirate the ethanol out and discard. Repeat for a total of two washes. Freshly made 70% ethanol may be used instead.

Note: It is important to perform these steps with the reaction plate on a 384 well magnetic plate. Do not disturb the separated magnetic beads. Be sure to remove all of the traces of ethanol from the bottom of the well as it is a known PCR inhibitor.

A drying time of 5 min at room temperature is optional to ensure all traces of ethanol are removed. Take care not to over dry the bead (beads appears cracked) as this will significantly decrease DNA yield.

6. Off the magnet plate, add 15 µL of elution buffer (Reagent grade water, Tris-Acetate pH 8.0, or TE) to each well and pipette mix 5 times.

A 15 µL elution volume will ensure the liquid level will be high enough to contact the magnetic beads. A greater volume of elution buffer can be used, but using less than 15 µL requires extra mixing (to ensure the liquid comes into contact with the beads) and may not fully elute the entire product. Elution is quite rapid and it is not necessary for the beads to go back into solution.

When setting up downstream reactions, pipette the DNA from the plate while it is situated on a 384 well magnetic plate. This will prevent bead carry over (however, beads will not inhibit thermal cycling reactions).



* All trademarks are property of their respective owners.

Disclaimer from ALINE Biosciences:

This product is for research use only and has not been validated for diagnostic purposes.

ALL INFORMATION PROVIDED IN THIS DOCUMENT BY ALINE IS ON AN “AS IS” BASIS ONLY. ALINE PROVIDES NO REPRESENTATIONS AND WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTIES OF FITNESS FOR A PARTICULAR PURPOSE, MERCHANTABILITY AND NONINFRINGEMENT. IN NO EVENT WILL ALINE BE LIABLE TO ANY PARTY FOR ANY DIRECT, INDIRECT, SPECIAL OR OTHER CONSEQUENTIAL DAMAGES FOR ANY USE OF THIS DOCUMENT INCLUDING, WITHOUT LIMITATION, FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT OR OTHERWISE, EVEN IF WE ARE EXPRESSLY ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.