


Sequencing Dye Terminator Removal ~ Superior Results & Cost Effective

PureSEQ™ MP provides superior purification efficiency by utilizing proprietary paramagnetic bead technology to effectively remove dye terminators, primer-dimers, salts, and nucleotides for traditional Sanger Sequencing workflows.

Advantages

- Long Reads with High QV Scores
- Eliminates GC n+1 peaks
- BigDye usage can be as low as 1/32 of ABI recommendation
- Validated and implemented in Clinical Reference Labs
- Lower repeat rate

Comparison Data



Ordering Information

Catalog Number	Product Size
N-4001-5	PureSEQ MP, 5mL. 600~1,000 rxns
N-4001-50	PureSEQ MP, 50mL. 5,200~ 10,000 rxns
N-4001-250	PureSEQ MP, 250mL. 25,000~50,000 rxns
N-4001-500	PureSEQ MP, 500mL 50,000~100,000 rxns
N-4001-BULK	Largest size / Custom Bulk



Aline Biosciences
200 Innovative Way Ste. 1305
Nashua, NH 03062 USA
order@alineliosciences.com